## DE Shiatsu-Massagekissen MC-81E

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## GB

## Safety instructions

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

### Power supply

- Before connecting the device to your power supply, make sure that it is switched off and that the electrical voltage specified on the device's type label matches that of your mains power.
- Always make sure that the plug is easily accessible when the device is connected.
- Keep the the cable and the massage seat away from heat, hot surfaces, moisture and liquids. Never get hold of the power plug with wet hands.
- Do not take hold of a device if it falls into water. Pull out the power plug from the socket immediately.
- Always unplug the equipment immediately after use.
- To disconnect the unit from the power supply, always pull the plug out of the socket. Do not pull on the cable!
- Do not carry, pull or twist the massage seat on the cable.
- You must not continue to use the unit if the lead is damaged. For safety reasons, it may only be replaced by authorized service personnel. In order to prevent hazards, always send the unit to the service centre.
- Ensure that the cable will not cause tripping. The cable must not be kinked, trapped or twisted.

### For special people

- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Do not use this unit to supplement or replace medical treatment. Chronic complaints and symptoms may be made worse.

- You should refrain from using the Shiatsu massage cushion, or should consult your doctor before using it, if
  - you are pregnant,
  - you have a cardiac pacemaker, artificial joints or electronic implants,
  - you suffer from one or more of the following illnesses and complaints: Circulation disorders, varicose veins, open wounds, bruises, broken skin, inflammation of the veins.
- Do not use the appliance near your eyes or near any other sensitive areas.

- Take care when using the heater. Do not use if you have circulatory problems or on areas not sensitive to heat.

### Risk of burning!

- The surface of the appliance is hot. People who are insensitive to heat must be careful when using the appliance.
- If you experience any pain or discomfort during the massage, stop immediately and consult your doctor.
- Consult your doctor if you have any questions concerning the therapeutic use of massage cushions.
- Please ensure that the polythene packing is kept away from the reach of children! Risk of suffocation!

### How does the ecomed Shiatsu massage cushion work?

- Your MC-81E Shiatsu massage cushion is a state-of-the-art device for genuine shiatsu-neck and back massages. It is also suitable for massaging your feet and legs. The MC-81E Shiatsu massage cushion has four rotating heads **1** which rotate in opposite pairs. While the left pair, arranged on a rotating plate, turns clockwise, the right pair is also arranged on a rotating plate, turn anticlockwise. The combination of intensive shiatsu massage and radiated heat makes the massage cushion an effective and relaxing treatment.

### Operating the unit

- Place the shiatsu massage cushion **1** on the back of a high-backed chair and fasten it to the back with the strap.
- Sit on the chair and check that the cushion is in the position you want.
- Connect the power unit **2** and the massage cushion and plug in the power unit to the mains. Ensure that the power supply is easily accessible.
- Press the massage button **3** to activate/deactivate the massage function. When the unit is activated the massage heads **1** are illuminated.

- When you have finished using the device, switch it off and pull the plug out of the power socket.
- Cleaning and maintenance**
- Before cleaning the unit, make sure it is switched off and disconnected from the mains. Clean the unit once it has cooled down using a soft, damp cloth.
  - Under no circumstances should you use aggressive detergents, strong brushes, solvents or alcohol, as they may damage the surface of the unit.
  - When cleaning, never immerse the unit in water and make sure that no water gets into the unit.
  - Wait until the device is completely dry before using it again.
  - Under no circumstances should you use aggressive detergents, strong brushes, solvents or alcohol, as they may damage the surface of the unit.
  - You should ideally store the device in its original packaging and keep it in a clean, dry place.
- Disposal**
- This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your local authority or your supplier for information about disposal.

**ecomed\*** is a trademark of MEDISANA AG  
Jägerbergstraße 19  
41468 NEUSS  
GERMANY  
E-Mail: info@medisana.de  
Internet: www.medisana.de

Murray's Health & Beauty  
School Lane  
Chandlers Ford  
Hampshire  
UK SO33 4YN  
Tel.: +44/2380460600  
eMail: cservices@paulmurrayplc.co.uk

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