
medisana®

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
GB Instruction manual
Massage Chair RS 100 / RS 150 / RS 200

Device and controls

❶ On/Off button with ❷ LED indicator
❸ “M” Button for switching the operating mode, with ❹ LED indicator
❺ HEAT button with integrated LED indicator
❻ Detachable head rest
❼ Sofa
❽ Shiatsu massage cushion
❾ Vibration massage cushion
❼ USB ports, for charging the massage cushions

Explanation of symbols

Do not puncture any of the components!

Not suitable for children under 3 years!

The device must only be used indoors!

Do not wash!

Do not bleach!

Do not dry the article in a tumble drier!

Do not iron!

Do not dry clean!

IMPORTANT
Non-observance of these instructions can result in serious injury or damage to the device.

WARNING
These warning notes must be observed to prevent any injury to the user.

CAUTION
These notes must be observed to prevent any damage to the device.

NOTE
These notes give you useful additional information on the installation or operation.

LOT number

Manufacturer

Recycling symbols / codes:
These are used to provide information about the material and its proper use and recycling.

CE

GB Important Instructions Retain for future use

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

- The device must only be used for its intended purpose as described in the instruction manual. Using the device for any other purpose invalidates the warranty.

Intended use

The device is used to massage certain body parts for relaxation and well being. It is not intended for medical therapeutic use. It is not suitable for clinics, massage practices, wellness studios etc. The massage device is not a medical device and is not intended for the treatment of pain. In case of doubt or uncertainty, consult your doctor before use.
The item consists of the actual sofa ❼, the Shiatsu massage cushion ❽ and a vibration massage cushion ❾.
The sofa ❼ is without any massage functionality.
The Shiatsu massage cushion ❽ and vibration massage cushion ❾ are working independent of each other, they have the same type of control panel ❶❷❸❹❺, and each is supplied with a rechargeable battery and a USB charging port. Therefore we will simply speak of “massage cushion”, “device” or similar in this instruction manual, and unless explicitly specified, the descriptions apply to ❽ and ❾ in the same way.

- **The device is not to be used for commercial use or medical purposes.**
- **If you have health concerns, consult your doctor before using the device.**

Only use the device indoors.

Do not use the device in damp rooms (e.g. when bathing or showering).

Safety Instructions – About the power supply

- Do not carry, hold or drag the device by the cable!
- Batteries must never be squeezed, disassembled or affected by hot temperatures (above 60°C). There is the danger of fire, explosions and the risk of a leakage of dangerous and poisonous chemicals.
- This device contains non-replaceable batteries.
- The batteries contained in this device may be replaced only by service technicians that have been authorised by **medisana**.
- This device must only be supplied at safety extra low voltage (SELV) corresponding to the marking on the appliance.
- The battery life gets reduced if it is not used for a long time. It is therefore recommended to charge the battery at least once or twice per month.
- Use a USB mains adapter with an output of 5 V and (at least) 2 A for charging the battery. A larger output voltage will damage the device.
- Before connecting the USB mains adapter to the power supply, please ensure that the input voltage stated on the rating plate is compatible with your mains supply.
- Only insert the power supply unit into a socket when the device is switched off. **The device must not be used while it is being charged.**

Safety Instructions – Specific groups of people

- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children must be supervised to ensure that they do not play with the device.
- Do not use this unit to supplement or replace medical treatment. Chronic complaints and symptoms may be made worse.
- Do not use the device, or only after consulting your doctor, if:
 - »you are pregnant,
 - »you have a cardiac pacemaker, artificial joints or electronic implants,
 - »you suffer from one or more of the following illnesses and complaints: Circulation disorders, varicose veins, open wounds, bruises, broken skin, inflammation of the veins.
- If you experience any pain or discomfort during the massage, stop immediately and consult your doctor.
- The surface of the appliance is hot. People who are insensitive to heat must be careful when using the appliance.
- Before using the massager, consult your doctor regarding any unexplained pain, if you are undergoing medical treatment or using medical equipment.
- Children under the age of three must not use this device due to their inability to react to overheating.

Safety Instructions – Using the device

- Always carefully inspect the device for damage before use. Do not attempt to use a defective device.
- Do not use the device if it shows signs of damage, if it is not in perfect working order, or if it has been dropped or become damp. In order to prevent hazards, always send the unit in to the service centre.
- Do not use the device before going to sleep. The massage has a stimulating effect.
- Never use the device while driving a vehicle or operating a machine.
- The maximum operating time for one session is 15 minutes. Prolonged use reduces the service life of the device.
- Very long periods of use can lead to overheating. Always allow the device to cool down before using it again.
- Never put the unit down next to an electric stove or other sources of heat.
- Do not allow the massage cushion to come into contact with pointed or sharp objects.
- Do not cover up the appliance when it is switched on. Never use it under blankets or cushions. There is a risk of fire, electric shock or injury.

Safety Instructions – Maintenance and cleaning

- The device is maintenance free. Do not try to open it. There are no parts inside that are to be cleaned or serviced by the user.
- You may only clean the device yourself. If the device malfunctions please do not repair it yourself as this leads to the automatic expiry of all warranty rights. Consult your specialist dealer and have repairs carried out by the service centre to avoid dangers.
- Do not immerse the device in water or other fluids. If, nonetheless, liquid should penetrate the device, unplug the mains cable immediately.
- Cleaning and user maintenance must not be carried out by children without supervision.

Scope of delivery

Please check first whether the unit is complete and does not show any damage.
If in doubt, do not operate the unit and contact your dealer or service centre.
The scope of delivery includes:
- 1 **medisana** Massage chair **RS 100 / RS 150 / RS 200**, consisting of:
 - »Sofa, detachable head rest,
 - »Shiatsu massage cushion, vibration massage cushion.
- 1 Charging cable, USB-C to USB-C(cable length: approx. 1 m)
- 1 Quick Start Guide (short instructions)

If you notice any transport damage when unpacking, please contact your dealer immediately.

Packaging is recyclable or can be returned to the raw material cycle. Please dispose of packaging material that is no longer required properly.

WARNING
Please ensure that the packaging material is kept away from the reach of children! Risk of suffocation!

Charging the battery

- Use a 5V mains adapter with an output current of (at least) 2 A. Connect the supplied USB charging cable to the adapter and to the USB port ❿ of the device.
- The device must not be used during the charging process.
- The LED ❷ above the On/Off button indicates the charging process by flashing in red colour. When the battery is fully charged, the LED ❷ is lit continuously in green colour.
- If the LED ❷ is lit in red colour during operation, then the remaining charge is low – the battery should be recharged as soon as possible. When the remaining charge is almost empty the LED ❷ starts flashing in red colour, and the device switches off shortly afterwards.

Using the device

- Press the On/Off button ❶ for (approx.) 2 seconds to switch the device on. It starts immediately with the massage. The operation is indicated by the LED ❷ being lit in green colour. (If the battery charge is getting low the colour changes to red, see previous section.)
- Press the Heat button ❺ for (approx.) 2 seconds to switch the heat function on. The button ❺ is lit in red colour when the heat function is active. Press the Heat button ❺ again for (approx.) 2 seconds to switch the heat function off.
- A short press on the “M” button ❸ changes between two operating modes (which is indicated by the LED ❹ flashing twice in blue colour):
 - »Shiatsu massage cushion ❽: The “M” button ❶ changes the massage rotation. After switching on, the “forward” direction is active by default.
 - »Vibration massage cushion ❾: The “M” button ❶ toggles between wave vibration and alternating left/right vibration. After switching on, the wave vibration is active by default.
- Press the On/Off button ❶ again for (approx.) 2 seconds to switch the device off. If you do not switch it off manually it switches off automatically after (approx.) 15 minutes.
 - »Shiatsu massage cushion ❽: With a fully charged battery, the device can be used for approx. 30 minutes, depending on whether heat is used.
 - »Vibration massage cushion ❾: With a fully charged battery, the device can be used for approx. 80 minutes, depending on whether heat is used.

Cleaning and maintenance

- Before cleaning the unit, make sure it is switched off and disconnected from the mains.Remove the charging cable.
- Before cleaning the appliance, make sure that the appliance is switched off and the power supply unit is disconnected. Allow the unit to cool down.
- Clean the device with a slightly damp sponge only. Never use brushes, aggressive cleaning agents, petrol, thinners or alcohol. Rub the unit dry with a clean dry cloth.
- May not be bleached, ironed or dry cleaned!
- Never immerse the appliance in water to clean it and ensure water cannot get into the appliance.
- Wait until the device is completely dry before using it again.
- Finally, reinsert the device into its correct place within the sofa.

Technical Data

Name and model:	**medisana** Massage chair **RS 100** / **RS 150** / **RS 200**
Power supply (Shiatsu massage cushion ❽):	Li-Ion battery 7.4V, 2500mAh; input: USB-C, 5V⎓ 2,0A; charging time: approx. 3 hours
Power supply (vibration massage cushion ❾):	Li-Ion battery 7.4V, 2500mAh; input: USB-C, 5V⎓ 2,0A; charging time: approx. 3 hours
Charging time:	approx. 3 hours (each for ❽, ❾).
Max. operating time:	15 minutes, with automatic switch-off
Maximum carrying capacity:	150 kg
Dimensions:	approx. 99 × 98 × 67 cm (without head rest ❻)
Weight:	approx. 22 kg
Operating conditions:	Only use in dry rooms at −5°C to +40°C
Storage conditions:	In a cool, dry place

Name	Colour	Item number / EAN number:
RS 100	light grey	87751 / 40 15588 87751 3
RS 150	dark grey	87752 / 40 15588 87752 0
RS 200	orange	87753 / 40 15588 87753 7

Disposal

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

medisana GmbH,
Carl-Schurz-Str. 2, 41460 NEUSS, GERMANY

In accordance with our policy of continual product improvement, we reserve the right to make technical and visual changes without notice.

The current version of this instruction manual can be found under https://docs.medisana.com/87751.
Information about service can be found here:
www.medisana.com/servicepartners

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